

Final Design Report
Elbow Motion Simulator

December 6th, 2024

The University of Western Ontario

Department of Mechanical and Materials Engineering

MME 2259A – “Product Design and Development”

Dr. Louis Ferriera

Final Project Documentation

Executive Summary

Our project addresses the need for total elbow replacement implants due to the increase in elbow replacement surgeries. We did this by designing an elbow implant motion simulator. This simulator is designed to replicate the biomechanical functions of a human elbow and provide a way of testing elbow implants in order to determine their reliability for medical device manufacturers and healthcare professionals.

The simulator uses a cadaveric specimen to assess the performance of an elbow implant. The design focuses on precision, durability and ease of use to ensure that it meets customer needs. Features include actuators mimicking the motion, the load bearing capacity of up to 1kN and chrome stainless steel for corrosion resistance and strength to help with its long-term durability.

Our specification and planning phase included using a Gantt chart, Critical path method, and Quality Function Deployment to ensure efficient use of our time in creating a design that meets the needs of the customer. The concept generation phase allowed us to brainstorm ideas to decide the best way to simulate the natural functions of the elbow. The final design was selected through the Go/No-Go matrix to ensure we picked the best design that met both the customer's needs and functional requirements.

The simulator was designed to create a standardized way of testing elbow implants to ensure patients are receiving the best possible plant that meets their needs. This design ultimately is intended to improve the comfort, functionality and safety of an implant and meet the requirements of a typical medical device.

Introduction

There is an increased number of elbow replacement surgeries, and there are no widely marketed elbow implant testing simulators. As the elbow implant demand increases it is necessary to set a standard of the range of motion, durability and how long each implant would last. The primary purpose of the elbow simulator is to evaluate the biomechanical performance of elbow implants on cadaveric specimens. Testing an elbow implant is necessary to determine its durability, safety and effectiveness. This simulator would benefit medical device manufacturers and medical professionals as it would provide them with a reliable machine for testing elbow implants, ensuring the safety and comfortability of a patient.

Project Plan

Phase 1: Specification and Planning

In this phase we established the problem statement, by identifying gaps in the market. We outlined our constraints to ensure the product effectively meets the needs of the customer. Planning the tasks for the project was important as it ensured each step is built upon the previous one, for the project to be completed efficiently.

Customer Needs

- **Realistic Joint Mechanics** - literature indicates that the elbow's natural range of motion includes flexion, which allows movement between 130 and 154 degrees, and extension, which ranges from -6 to 11 degrees. Additionally, the elbow facilitates pronation, with a range of 75 to 85 degrees, and supination, spanning 80 to 104 degrees so in our elbow simulator we will strive to get the motion as close as possible to these values. [1]
- **Precision and accuracy in motion** - Actuators will deliver smooth and controlled movements that closely mimic human biomechanics
- **Load Testing** - literature shows that a joint load in an elbow implant can induce peak values of about 1kN so our simulator should be capable of withstanding and testing that amount of load. [2]
- **Flexibility in configuration** - the design should have adjustable components that allow for customization by the client however they see fit
- **Ease of Use** - The design must be easy to set up and operate while reducing the time and effort required
- **Safety and regulatory compliance**
- **Long Lasting:**
- **Dynamics Load Response:** The design must be able to accurately detect and show load forces while simulating the natural movement of a forearm
- **Durability:** The simulator must be able to withstand repeated use in realistic conditions (stress and movement). Material must also show little to no wear during and after each testing is done.
- **Compact/reasonable size:** The design must be sturdy and rigid, with little to no wobbling when the cadaveric arm flexes, tenses, supinates and pronates.

Problem Statement

With the increased need for elbow implants, it is important to have standardized testing and performance benchmarks. When patients opt for a total elbow prosthesis, they expect it to restore their original range of motion that allows them to live with little to no pain. Correctly evaluating implant designs ensures that they meet the necessary biomechanical functions, to provide patients with reliable and effective results. Based on this information we set our problem statement as: **Current ex-vivo elbow simulators lack precision and need accuracy to replicate the mechanics and movements of a human elbow.**

Product Design Specification (PDS)

The PDS assisted us in determining the things we need to consider before concept generation, to ensure it meets our customer needs. Extensive research needed to be done to determine the necessary biomechanical functions and safety standards, which is critical for accurately testing a medical device.

Specification Type	Performance	Operating Environment	Standards	Materials	
	The elbow simulator must perform the natural functions of a human elbow.	Elbow simulator must withstand pressure ranges of at least 6MPa as contact pressures on cartilage can reach up to 5.5MPa.	ASTM F ASTM E8/E8M ISO 13485 ISO 14708	Not yet determined but should be lightweight and have corrosion resistance and must be able to withstand heavy weight from implant testing procedures.	
		In an environment with atmospheric pressure of 101.325 KPa.			
Specification Type	Ergonomics	Aesthetics, Appearance and Finish	Competition Benchmarking	Quality and Reliability	Testing and Inspection
	Easy to setup and mount the arm with the necessary wires	Simulator must look clean and polished; any wires must not interfere with the view of part inserted in the simulator	Competition can only test implants that slide against each other (i.e. joint implants)	High durability so it can last for thousands of cycles	Regular sanitization of apparatus to get rid of residue from previous cadaveric specimen
			Don't account for the other muscles that act on the joint	Be compatibility with various elbow implants.	Ensure that the simulator is able to handle certain amounts of weight after multiple uses, by testing the machine

[illegible]

Specification Type	Manufacturing facilities	Patents	Design Schedule	Company Constraints	Social and Political Factors	Safety
	Must be produced in a facility that means the medical device manufacturing standards	Joint simulator testing machine	About 20 days for part modelling and Assembly of the simulator to ensure movement of an elbow is correctly mimicked	The product is niche and maybe be outperformed by competitors different joint implants	Must improve the quality of implants	Wiring should be properly insulated
	Be in a facility that is near where we source some or all the materials to reduce sourcing and shipping cost	Smart Joint Implant Sensors			Must comply with sustainability and environmental regulations and laws	Most have disposable or replaceable components for parts that interact with the cadaveric specimen

Figure 1- Product Design Specification

Work Breakdown Structure

The WBS assisted in dividing the project into manageable components. It was important to divide the tasks into specification and planning, concept generation and the design phase to show how each stage builds on to the next ensuring an efficient workflow with designated responsibilities for each team member. This was used to assist in the preparation of the critical path method and Gantt chart, to ensure every task was accounted for.




Figure 2 – Work Breakdown Structure

The Quality Function Development (QFD) was utilized to translate the identified customer requirements and how it will relate to our technical requirements. This model helped us see how certain solutions or ways we would satisfy those customer needs related to each other and whether certain technical requirements that were needed to satisfy customer needs were doable or out of reach. Additionally, the QFD helped us compare advantages and disadvantages we had over the competition and see how well our design would do against them.



Gantt Chart

Gantt Chart was used to set deadlines and gave us an indication if we were on track during each task. The durations gave us a basis on how long we should take on each task. The chart was designed in a way to give us flexibility in case of delays. The schedule allowed us to make the entire process more efficient and ensured each team member knew what to do next.




Figure 4 – Gantt Chart

Critical Path Method

This was used in tandem with the Gantt to ensure we are on the right track and let us determine the maximum and minimum project duration. This allowed us to plan the best course of action, to ensure we had ample time to finish the long components of our project, specifically the design phase. The CPM allowed us to prioritize which task, to ensure we complete everything by the deadline.




Figure 5- Critical Path Method

Concept Generation

During the initial stages of ideation, we begin with a single concept inspired by our analysis of our competitors. We adapted and improved the original concept which enabled us to brainstorm and design 8 ideas. This strategy allowed us to refine the original idea, creating a design with increased functionality, so we could differentiate ourselves from the competition.

Sketches

To better understand and evaluate our ideas, we created sketches for each concept (see Figures 8-15). These visual representations provided a clearer understanding of how each concept would function and helped identify potential challenges early in the process. Even concepts that we initially believed were not feasible were sketched, to allow us to compare and analyze all the ideas effectively. These sketches played a crucial role in our decision-making process, enabling us to narrow down and select the most promising concept to move forward with.




Figure 6 – Guide rails to allow more flexion accuracy + pneumatic actuator connected to ulna and radius to simulate elbow rotation + Clamps to fix humerus in position




Figure 7 – Explanation of guide rail mechanism




Figure 8 – Add more alignment systems to better align the linear actuators with each other and their respective muscle




Figure 9 – Actuator pushes the arm up and down at the hand, which exerts a force on the elbow. Force sensors that are surgically added in the elbow measure this force




Figure 10 – Add more elastics to give a more realistic representation of the muscles in the elbow area




Figure 11 – A clamp to hold down the humerus bone, with movement through pneumatic cylinders




Figure 12 – Movement through pneumatic cylinders and forearm rotation through a rotary actuator around the radius and ulna bones.




Figure 13 – Flexion and extension through a pulley system with weights at one end and forearm rotation through a rotary actuator around the radius and ulna bone

Functional Decomposition and Morphological Analysis

During our concept generation stage, we employed functional decomposition and morphological analysis to structure our design process. As time went on, our designs became more creative and complex. We focused on the main functions we wanted our simulator to perform and brainstormed methods to achieve them for each concept. For example, in one of our concepts, we wanted to incorporate supination and

pronation in the forearm. To do this, we thought to use a rotary actuator to rotate the ulna and radius, the bones in the forearm, at the elbow joint, which is how supination and pronation occur in the human body. The combination of these two design tools allowed us to generate ideas with a more feasible implementation.

Concept Selection – Go/No-Go Screening

To decide on the concept, we used a method called the Go/No-Go Screening. This method compares each concept to a metric, which in our case was the customer needs. The comparisons that met the criteria received a 1, which meant ‘go,’ the ones that may meet the criteria received a ‘0,’ which meant ‘maybe,’ and the ones that did not meet the criteria received a –1, which meant ‘no-go.’ The concept that had the greatest sum, that is, had the greatest net positive comparisons, was the concept we decided to pursue (see Figure 15). (See Appendix A for the final design)

	Precision and Accuracy in Motion	Load Testing	Flexibility of Configuration	Ease of use	Safety	Long Lasting	Dynamics Load Response	Durability	Realistic Joint Mechanics	Compact/ reasonable size	Σ
Clamp + mount to hold down humerus bone + movement of muscles through pneumatic system	0	-1	-1	1	1	0	1	1	-1	1	2
Guide rails to allow more flexion accuracy + pneumatic actuator connected in between the ulna and radius through a rope to flex and extend the forearm, simulating the elbow joint movement+ clamps to fix the humerus bone in place	1	1	0	1	1	0	1	1	0	1	7
Add more elastics to give a more realistic representation of the muscles in the elbow area	1	1	-1	1	-1	-1	1	0	0	1	2
Add more alignment systems to better align the linear actuators with each other and their respective muscle	1	1	-1	1	1	1	1	0	0	1	6
Actuator will connect to the tendon or cable attached to the biceps tendon. when the actuator contract or expand, it pulls on the tendon which simulate the contraction and tension of the biceps muscle (moving the forearm up and down). + Connection to ulna and radius of rotary actuator to simulate elbow rotation	1	0	-1	1	1	0	1	0	1	1	5
Actuator pushes arm up and down at the hand, which exerts a force in the elbow. Force sensors that are surgically added in the elbow measure this force	-1	1	-1	1	1	0	1	-1	-1	1	1
Pulley system where thin and strong cable is connected to forearm and counterweight is at the other end of the cable, moving the forearm up and down. Additionally it has an electric rotary actuator connected in the ulna and radius	-1	1	1	1	0	0	1	1	0	-1	3

Figure 14 – Go/No-Go Screening

Design Review

Frame Support

This is attached to the frame post as a way for the simulator to have structural stability, by bearing the load of the entire mechanism. The stability it provides allows the user to maintain the safety necessary for using the simulator.

Material: Chrome Stainless Steel

Used for its corrosion-resistant properties, as we mentioned in the PDS. With lubrication and constant sanitation, this material will increase its durability and lifespan.





Figure 15 – Frame support CAD and CAD drawing

Frame Post

This is an important part as it assists in the ergonomics of our design as it allows for easy setup of the simulator as it is used to place the rotary actuator, linear actuator, and the cadaveric specimen. The frame post does not have a backboard as it allows for the motion to be clear and visible.

Material: Chrome Stainless Steel

It is used for its corrosion-resistant properties, as well as its low electrical conductivity, therefore increasing the safety if it comes in contact with electrical components.




Figure 16 – Frame Post CAD drawing and CAD Drawing

Hinge Mount:

This part is connected to the sheet metal and with the *door hinge* part. The hinge mount, as the name suggests, acts as a mount for the door hinge that will be placed inside this part.

Material: Chrome Stainless Steel is used for its high yield strength and its resistance to corrosion and wear so that it can hold the door hinge part for thousands of cycles and still work perfectly fine.




Figure 17- Hinge CAD and CAD Drawing

Door Hinge:

The door hinge allows for an opening of 90 degrees in the clamp to allow for the humerus bone to be placed and then clamped down.

Material: Chrome Stainless Steel is used for its high yield strength and its resistance to corrosion and wear so that it can hold the door hinge part for thousands of cycles and still work perfectly fine.




Figure 18- Door Hinge CAD and CAD Drawing

Frame Plates:

This is used to connect the frame supports and frame posts, it assists in increasing the stability of the simulator and therefore its reliability. Which is necessary for the heavy weight of the other components such as the actuators and cadaveric specimen.

Material: AISI 1020 Steel is used for its high corrosion resistance and its high strength.




Figure 19 – Frame Plates CAD and CAD Drawing

Cadaver Holder

This is used to hold up the bone and linear actuator, so the bone is visible and clear when performing the necessary testing procedures. It allows for the bone to move without it being moved to unwanted positions.

Materials: Chrome Stainless Steel as it has high strength, to be able to hold the weight of both the entire cadaver and linear actuator.




Figure 20 – Cadaver Holder CAD and CAD Drawing

Forearm Dynamics Stabilizers

This is used to stabilize the elbow, so it can flex and extend without the elbow moving side to side. This ensures that during testing procedures, flexion and extension can be accurately measured.

Material: Natural Rubber is used for its elasticity, as it mimics the function of the tendon. So it allows for the controlled movement of the elbow during testing procedures.




Figure 21- Forearm Dynamic Stabilizers CAD and CAD Drawing

Conclusion

The development of the elbow motion simulator is necessary for advancing the testing procedures of the elbow implants. By assessing all the product design specifications, we were able to determine the necessary precision and durability while maintaining safety, to replicate elbow mechanics of flexion and extension. This simulator sets benchmarks for the testing of implants by inspecting their range of motion, load bearing capacity, and longevity of implants. The simulator ensures that the implant complies with the safety standards of medical devices.

We were guided by the customer needs and used specification and planning methods such as a WBS, Gantt chart and CPM to lead us to comprehensive design. We exhausted design ideas to design a simulator that meets both the customer needs and medical device standards.

Each design component, from the frame support to the rotary actuator mount, was designed to meet the structural needs of the simulator to ensure the stability and reliability of the design to extend its lifespan. The material used also assisted in the reliability of the design by using chrome stainless steel and natural rubber, so the biomechanical functions of the elbow can be replicated with as few problems as possible.

In conclusion, the elbow implant motion simulator is designed to become an essential tool in testing implants. It will contribute to developing safer and more reliable implants, ensuring patients can better their quality of life with better use of their elbow.

Recommendations

- An actuator to simulate the triceps could have been more realistic when simulating the flexion of the forearm, as the current design replaces the flexion done by the triceps with the force of gravity.
- An extra shaft or adding an extra smaller bicep actuator to simulate the tendon that would simulate the supination movement done by the bicep.

-Adding a replacement for the hand at the end of the forearm so that it can simulate the elbow joint moving when heavily loaded

-Many of the parts are welded which doesn't allow for disassembly


- There is only one clamp on the humerus, which may not be very stable when the arm is heavily loaded. More clamps or ways to hold the humerus in place should be added.

Appendix

[A]


Front view

frame



Back View

frame



References

- [1] E. L. Zwerus et al., “Normative values and affecting factors for the elbow range of motion,” Shoulder & elbow,
<https://pmc.ncbi.nlm.nih.gov/articles/PMC6555111/#:~:text=The%20elbow%20joint%20allows%20us,80%C2%B0%20to%20104%C2%B0>. (accessed Sep. 18, 2024).
- [2] S. K. Armah, “Stress analysis of an artificial human elbow joint,” Stress Analysis of an Artificial Human Elbow Joint: Application of Finite Element Analysis,
<https://thescipub.com/pdf/ajeassp.2018.1.18.pdf> (accessed Sep. 18, 2024).
- [3] Combs, Taylor et al, “Active Motion Laboratory Test Apparatus for Evaluation of Total Elbow Prostheses,” Journal of Hand Surgery Global Online, Volume 6, Issue 1, 21 – 26
[https://www.jhsgo.org/article/S2589-5141\(23\)00143-3/fulltext#fig3](https://www.jhsgo.org/article/S2589-5141(23)00143-3/fulltext#fig3) (accessed Nov.14, 2024).
- [4] Lazar, Michael & Lurski, Alex & Ghafurian, Soheil & Chen, Linda & Uko, Linda & Tan, Virak & Li, Kang. (2015). In vitro Dynamic Simulation of Elbow Motion. 10.1109/NEBEC.2015.7117109.
https://www.researchgate.net/publication/278849845_In_vitro_Dynamic_Simulation_of_Elbow_Motion (accessed Nov. 14, 2024)